

mobiFUSER

Home care infusion pump



**Elastomeric
Infusion Device**

*Instructions
For Use*



Elastomeric Infusion device

This brochure is intended as a guide for the proper use of the **mobiFUSER** disposable infusion pump.

The **mobiFUSER** is a safe, effective and durable delivery device. The mobiFUSER delivers medication on a continuous basis over a specified period of time.

Please consult your treating health professional if you have any questions after reading this guide.



Introduction





The **mobiFUSER**



1. LANYARD

The clip can be used as an alternative to the carry case supplied with the **mobiFUSER**.



2. OUTER CASE

An outer cover to protect the medication reservoir.



3. RESERVOIR

The reservoir or balloon holds your medication. The slow deflation of the balloon pushes the medication through the tubing to your catheter.



4. INFUSION VOLUME INDICATOR

The black line at the top of the reservoir indicates remaining volume.



5. ACCESS PORT

A port where the medication is added.
DO NOT ACCESS THIS PORT.




Parts



6. ADMINISTRATION TUBING

Delivers medication from the balloon to your catheter.



7. AIR FILTER

This filter traps any air in the tubing.




8. FLOW RATE CONTROL

The coloured tubing restricts flow of medication to a specific rate.



9. CONNECTOR/CAP

The connector at the end of the tubing that attaches to your catheter. Comes capped. Keep this cap to seal your used **mobiFUSER** when the infusion is finished.



mobiFUSER Storage

The **mobiFUSER** is usually stored in a refrigerator (2 - 8°C) or at room temperature and (where possible) out of direct light. Store the product as per label instructions. Or consult your treating health professional for storage instructions.

IF STORED IN A REFRIGERATOR

- Store **mobiFUSERs** in a clean separate compartment away from food
- Store **mobiFUSERs** in the protective outer pouch provided
- Allow enough room around the **mobiFUSERs** for cold air circulation
- Prior to use, remove **mobiFUSER** from refrigerator and allow to warm to room temperature, protected from light (Maximum 4 hours prior to use, consult your treating health professional for the most appropriate time frame)

IF STORED AT ROOM TEMPERATURE

- Store in a clean environment away from light, moisture and heat
- Store in a safe location out of the reach of children and pets



Storage

Other Important Information To Know

ATTENTION

- The **mobiFUSER** is single use only, discard it after use as instructed by your treating health professional
- Make sure there are no bends or twists in the administration tubing line as this will reduce the amount of medication flowing through the line

CONSULT TREATING HEALTH PROFESSIONAL IF

- The **mobiFUSER** is damaged or leaking [keep the damaged device to show your treating health professional]
- Your catheter becomes dislodged
- The **mobiFUSER** has not infused medication over desired time period [flow rates can vary depending on your catheter]
- The fluid in the **mobiFUSER** has become discoloured or has any sediment
- Any information on the label is incorrect [e.g. your name, medication, dose, expiry date]



Changing Your **mobiFUSER**

CONNECTION

Always clean your hands before handling the mobiFUSER or your catheter.

1. Remove mobiFUSER from protective outer pouch.
2. Remove the connector cap and watch for fluid to slowly start forming a droplet at the end of the connector. Replace cap.
3. Access the catheter 'aseptically' (using clean technique) as instructed or as per local procedures. If required, test flush catheter with appropriate sterile fluid to check that it flushes well.
4. Remove cap and connect to catheter (turn connector clockwise).
5. Secure mobiFUSER to skin with tape placed in-between air filter and flow rate controlled tubing.
6. Place mobiFUSER in carry case provided.

Note: The connector cap can be saved to recap mobiFUSER when the infusion is finished.

Connection/Disconnection

DISCONNECTION

Always clean your hands before handling the mobiFUSER or your catheter.

1. Remove mobiFUSER from carry case.
2. Remove skin securing tape.
3. Disconnect connector by turning anticlockwise.
4. Access the catheter aseptically as instructed, or as per local procedures. If required, flush catheter with flush solution.
5. Recap used mobiFUSER with old connector cap.



Using the **mobiFUSER**

CHECKING INFUSION IS WORKING

The **mobiFUSER** delivers medication continuously over a desired period, flow rate may appear slow [i.e 1 drop every 15-20 seconds is standard]. The reservoir (balloon) gradually gets smaller over the infusion period. As the balloon deflates the black line [volume indicator] slowly moves downwards. You can check how much of your medication has infused by lining up the black line with the amounts on the outer case.

DAY-TO-DAY USAGE

SHOWERING/BATHING

The **mobiFUSER** should be kept dry at all times.

However, this should not stop you showering/bathing. Your treating health professional can instruct you on how to avoid getting the **mobiFUSER** wet when showering/bathing.

Usage



EVERYDAY LIVING

- Keep the **mobiFUSER** in the carry case provided, avoid exposure to prolonged periods in light
- Avoid any activity that may cause damage to the **mobiFUSER**
- Avoid exposure to temperature fluctuation.

If you have any questions about what you are able to do with the **mobiFUSER** connected, consult your treating health professional for suitable activities. For example, if you need to travel at altitude.

SLEEPING

- The **mobiFUSER** should be placed on top of your blankets/ quilt - NEVER under the covers
- Ensure the **mobiFUSER** has been securely attached to the catheter
- Always place the **mobiFUSER** close to the level of the catheter. Infusions rates may vary if the **mobiFUSER** is placed too high above or below the level of the catheter for prolonged periods [storage in the carry case is acceptable]





sladehealth.com.au